




XXXVI Reunión del grupo de electroquímica de la
Real Sociedad Española de Química
XVII Encontro Ibérico de Electroquímica

Research Project - Electrochemical science and technology



Screen printing of electrochemical biosensors on flexible substrates. Development of electronic prototype devices




Hugo Flávio Cunha da Silva

Supervisor : María Julia Arcos Martínez


Main Goal



- ❖ Development of screen-printed systems, using conductive textiles modified with enzymes and/or with nanomaterials, to be used as an autonomous biosensor without the need of a connected conventional potentiostat.




Wearable autonomous biosensor



Adapted from Malzahn K., *et al.* Analyst, 2011, 136, 2912



Screen-Printed Electrodes



Advantages:

- ✓ Automated, easy and low cost fabrication
- ✓ Reusable or disposable devices
- ✓ Possibility of electrode modification with bio and nanomaterials
- ✓ Miniaturized for fast analysis of small volume samples
- ✓ Portability for use in situ
- ✓ Reproducibility

Screen-Printed Electrodes



Flexible Substrates


- ❖ Incorporation into clothes,
- ❖ Incorporation into food packages,
- ❖ Noninvasive devices for “*in situ*”
- ❖ Incorporation of potentiostat



Screen-Printed Electrodes



Flexible Substrates



Inert and resistant properties

Mechanical resiliency, plasticity, thermal stability


Hydrophobicity, suitable dielectric and insulating properties

- Polyimide (Kapton),
- Polyethylene naphthalate (PEN),
- Polyethylene terephthalate (PET),
- Polytetrafluoroethylene (Teflon),
- Tattoo paper

- Gore-Tex
- Neoprene
- Polyester
- Nylon

Textiles

❖ L-lactic Acid



- Sports medicine

Training status and physical condition of the individual

- Clinical diagnosis

Respiratory insufficiency, shocks, heart failure and metabolic disorder

- Food industry

Indicator of fermentative process, storage quality and products freshness.

Analytical methods

- ✓ High-performance liquid chromatography
- ✓ Fluorometry
- ✓ Magnetic resonance spectroscopy
- ✓ Enzymatic Kits




Disadvantages

- × Time-consuming,
- × Require the sample pre-treatment
- × Expensive
- × Require educated personnel
- × Need of sample collection

Interest Analyte



❖ L-lactic Acid in sport medicine



Blood Levels

Rest - 0.5-1.5 mM

During Exercise - as high as 12 mM

Extreme conditions and extensive physical exercise - 25 mM

Interest Analyte



❖ L-lactic Acid in sport medicine

Samples	Blood	Interstitial fluid	Tears	Exhaled breath	Saliva	Sweat*
[LA]	0.5-20.0 mM	1.9-2.2 mM	1.0-5.0 mM	13.5-22.0 µM	0.1-2.5 mM	Max. 25.0 mM

L.Rassei et al. Anal Bioanal Chem, 2013, 406: 123-137; * J. Green et al. Eur J Appl Physiol (2004) 91: 1-6.

Sweat lactate

Sweat Eccrine gland energy metabolism



Exercise intensity



[Lactate]



Perspiration

Interest Analyte



❖ L-lactic Acid in sport medicine

Samples	Blood	Interstitial fluid	Tears	Exhaled breath	Saliva	Sweat*
[LA]	0.5-20.0 mM	1.9-2.2 mM	1.0-5.0 mM	13.5-22.0 μ M	0.1-2.5 mM	Max. 25.0 mM


L.Rassei et al. Anal Bioanal Chem, 2013, 406: 123-137; * J. Green et al. Eur J Appl Physiol (2004) 91: 1-6.




Interest Analyte





Lo_x - Lactate Oxidase




Electrode modification with LO_x provides selectivity in body fluids analysis.



+



+



Lactic Acid


Pyruvate

Hydrogen Peroxide

Summary Research Schedule


XXXVI

G1




Screen printing

G2





L-Lactic acid determination

- Nanomaterials
- Mediators




Conventional Substrates

G3



Flexible Substrates

G4



Potentiostat




- HPLC

Validation

Screen-Printing Electrodes




Screen-printing inks



Screen-Printing Electrodes



Flexible Substrates



Substrates

cellulose acetate

~~PET (18 μ m)~~ Adhesion

Gore-Tex

~~Nylon~~


~~Poliéster~~ Hydrophobicity

- ✓ Excellent inks adhesion,
- ✓ Hydrophobicity
- ✓ Endurance with screen-printing process,
- ✓ Suitable conductivity

Screen-Printing Electrodes




Flexible Substrates



Substrates

cellulose acetate

Gore-Tex




Cyclic Voltamperogram for Ferrocyanide 5 mM in KCl (100 mM). Scan Rate of 50 mV s^{-1} .

- ✓ Redox behavior similar to the one observed in conventional substrates (PET 500 μm).



XXXVI Reunión del grupo de electroquímica de la
Real Sociedad Española de Química
XVII Encontro Ibérico de Electroquímica

Research Project - Electrochemical science and technology



Thank You

Acknowledgements:

